



SIEMENS



# CLINITEK Advantus<sup>®</sup> Urine Chemistry Analyzer

PEP to Go

# Table of Contents

1. System Overview

2. Quality Control

3. Sample Processing

4. Maintenance



# System Overview

## Key Components and Features

1. Touch Screen  
Software interface, used to operate the analyzer
2. Printer  
Automatically prints completed results
3. Strip Loading Station  
User loads strip test for analysis
4. Incubation/Read Station  
Incubates and analyzes/scans the strip test
5. Waste Bin  
Stores used test strips
6. Barcode Reader  
Scans information into the system



# System Overview

## Key Components and Features

1. Memory Card Slot  
Update software
2. Keyboard/Barcode Connector  
Connect external keyboard and hand held barcode reader
3. Interface Connectors  
Connects the analyzer to a hospital information system
4. Power Supply  
Provides power to the analyzer
5. Power Switch  
Turns the analyzer ON/OFF



# Quality Control

## Preparing Quality Control

1. Add 12 mL distilled water or DI H<sub>2</sub>O into an appropriately labeled tube with Siemens Chek-Stix<sup>®</sup> control strips for urinalysis
  - Prepare a Positive and Negative control
2. Remove a Chek-Stix<sup>®</sup> control strip
3. Place the control test strip into the appropriate tube
4. Cap the tube and invert gently for 2 minutes
5. Allow the tube to stand for 30 minutes
6. Invert the tube 1 more time
7. Uncap the tube and discard the control strip per laboratory procedure
8. Repeat steps for the other control

## Performing Quality Control

1. Select **Menu**
2. Select **QC**
3. Enter **ID Number**
  - QC Lot ID can be entered using a keyboard or scan it from a barcoded label
  - Select **QC0001** and **A-Z** to manually enter information
  - Select **Enter** to return
4. Enter the Expiration Date
  - The date can be entered using a keyboard or scan it from a barcoded label
  - The year and month can be entered manually
  - Select **Enter** to save the date
5. Select **Enter** to continue
6. Completely immerse all reagent pads of the strip test into the quality control solution
  - Immediately remove the test strip
  - Run the edge against the side of the container when removing
  - Caution: Do not blot the edge of the strip. This could affect results.
7. Place test strip onto the strip loading station
8. Repeat for the other control

# Sample Processing

## Performing a Urinalysis Strip Test

1. Select the appropriate test strip
2. Check/change sequence number (if applicable)
3. Enter a Tech ID (if applicable)
4. Enter the Color and Clarity
  - Use the cycle keys on the Ready/Run screen
  - Default Color/Clarity values are YELLOW and CLEAR
5. Prepare the test strip for analysis
  - Completely immerse all of the reagent pads into the specimen
  - Immediately remove the test strip
  - Run the edge of the strip against the container when removing
  - Place the test strip onto the loading station with the reagent pads facing up
  - Place to the right of the push bar
  - Ensure the end of the strip is against the back wall

Note: DO NOT blot the urine strips – this may affect results.



# Stored Results

## Recalling Patient Results

1. Select **Menu** from the Ready/Run screen
2. Select **Memory**
3. Recall Results
  - **All patient results** – the number of results will display to the right, a list of all stored patient results will be listed
  - **Last batch of patient results** – patient results between the last pause in testing and latest test will be listed
  - **Search for stored results** – search for a specific patient ID
4. Patient results can be printed, transmitted, and deleted from this screen

## Recalling Quality Control Results

1. Select **Menu** from the Ready/Run screen
2. Select **Memory**
3. Select **All QC Results**
  - The number of results will display to the right, a list of all stored QC results will be listed
4. QC results can be printed, transmitted, and deleted from this screen

# Maintenance

## Clean the Exterior and Touch Screen

1. Turn the analyzer power OFF
2. Clean the exterior of the analyzer using a damp cloth and a mild detergent
3. Clean the touch screen, with a soft nonabrasive cloth dampened with a mild detergent and water
  - Do not use bleach to clean the display



## Remove and Clean the Push Bar

1. Turn the analyzer Off
2. Tilting the bar slightly upwards, pull the bar straight out
3. Clean the push bar with warm water and a mild detergent
4. Rinse thoroughly and dry with a paper towel or soft cloth
5. To reinstall the push bar refer to the *Maintenance: Reinstalling* tab in this document





# Maintenance

## Remove and Clean the Fixed Platform

1. Turn the analyzer Off
2. Remove the waste bin container and discard the used test strips into an appropriate container
3. Remove the fixed platform by pull the entire assembly towards you
4. Clean the platform with warm water and a mild detergent
  - Rinse thoroughly and dry with a paper towel or soft cloth
  - Allow the calibration bars to air dry
5. Inspect the calibration bars for scratches, marks, or discoloration
  - Discard and replace the platform with a new one if the bars can not be cleaned
6. To reinstall the platform refer to the *Maintenance: Reinstalling* tab in this document



# Maintenance - Cleaning


## Remove and Clean the Holddown Plate

1. Turn the analyzer Off
2. Press upwards on the tab at the back of the plate
3. Pull the other end from its retaining hole
4. Clean the holddown plate with warm water and a mild detergent
  - Wipe with a damp cloth
  - If the holddown plate is extremely dirty, soak in warm water and a mild detergent to loosen the dried residue
  - Rinse and dry with a paper towel or soft cloth
5. To reinstall the holddown plate refer to the *Maintenance: Reinstalling* tab in this document



## Remove and Clean the Moving Table

1. Turn the analyzer Off
2. Pull the moving table towards you
3. Clean the moving table with warm water and a mild detergent
  - Rinse thoroughly and dry with a paper towel or soft cloth
  - Use care when drying around the pins
6. To reinstall the moving table refer to the *Maintenance: Reinstalling* tab in this document



# Maintenance - Reinstalling

## Reinstalling the Moving Table

1. Hold the table with the small rectangular tab facing to the back
2. Align the 2 grooves on the bottom of the table with the edges of the platform on which the table rests
3. Gently push the table in until you hear the tab latch into the hold position
4. Check that the table is secure



## Reinstalling the Holddown Plate

1. Position the holddown plate with the arrow side facing up and the arrow pointing towards the back
2. Place the pin on the front of the holddown plate into the hold at the front of the fixed platform
3. Align the tab at the back of the holddown plate with the slot at the back of the platform
4. Snap the holddown plate into place. Listen for a loud click, indicating proper installation
5. Ensure the white calibration bars are visible



# Maintenance - Reinstalling

## Reinstalling the Fixed Platform

1. Align the 2 grooves on the bottom of the fixed platform with the arms extending forward from the analyzer
  - The flanges on the sides of the holddown plate align just outside the read area cover. The top edge of the platform aligns just under the cover
2. Gently push the platform in as far as possible
  - Push past the ridge to correctly position the platform
  - Caution: Do not force the platform. You may damage the moving table or fixed platform if incorrectly positioned



## Reinstalling the Push Bar

1. Hold the push bar at the indented end
2. With the end slightly upward, insert the peg on the other end of the bar into the hold in the push mechanism
3. Lower the push bar into place

## Reinstalling the Waste Bin Liner

1. Place a new liner into the waste bin



Please note that the learning material is for training purposes only!

For the proper use of the software or hardware, please always use the Operator Manual or Instructions for Use (hereinafter collectively "Operator Manual") issued by Siemens Healthineers. This material is to be used as training material only and shall by no means substitute the Operator Manual. Any material used in this training will not be updated on a regular basis and does not necessarily reflect the latest version of the software and hardware available at the time of the training.

The Operator's Manual shall be used as your main reference, in particular for relevant safety information like warnings and cautions.

Note: Some functions shown in this material are optional and might not be part of your system.

Certain products, product related claims or functionalities (hereinafter collectively "Functionality") may not (yet) be commercially available in your country. Due to regulatory requirements, the future availability of said Functionalities in any specific country is not guaranteed. Please contact your local Siemens Healthineers sales representative for the most current information.

The reproduction, transmission or distribution of this training or its contents is not permitted without express written authority. Offenders will be liable for damages.

All names and data of patients, parameters and configuration dependent designations are fictional and examples only.

All rights, including rights created by patent grant or registration of a utility model or design, are reserved.

Copyright © Siemens Healthcare GmbH 2020

Siemens Healthineers Headquarters  
Siemens Healthcare GmbH  
Henkestr. 127  
91052 Erlangen, Germany  
Telephone: +49 9131 84-0  
[siemens.com/healthineers](https://www.siemens.com/healthineers)