Dilutions

Concentration = Particles/Volume

or,
$$C = P/V$$

For example, if 0.1 ml of culture is used for a spread plate and 150 colonies are counted:

$$C = 150 \text{ CFU}/0.1 \text{ ml} = 1500 \text{ CFU/ml}$$

Where: C = the concentration of colony forming units (CFU) in the culture plated

P = the number of CFU on the plate

V =the volume plated

However, what if the aliquot plated was from a 1/10 dilution of the culture? It seems obvious that the concentration of CFU in the original culture would be 10X higher than the value calculated above, or 15,000 CFU/ml. How can this type of relationship be expressed mathematically, in order to simplify making dilutions and performing calculations later?

If a suspension is diluted, the number of particles doesn't change:

$$P_1 = P_2$$

but, since the volume increases, the concentration (P/V) decreases. Based on the equation for concentration above:

$$P = CV$$

and,

$$C_1V_1 = C_2V_2$$

That is, since the number of particles is the same in each sample, the concentration in each sample (expressed, for example as CFU/ml) times the volume (ml) would be the same.

To calculate dilutions:

$$Dilution = \frac{Volume \ added}{Total \ volume} = \frac{Volume \ added}{Volume \ of \ blank + Volume \ added}$$

For example, if a culture contains 1 x 10⁹ cells/ml and 0.1 ml is added to 9.9 ml of a blank,

$$Dilution = 0.1 \text{ ml}/(9.9 \text{ ml} + 0.1 \text{ ml}) = 0.1 \text{ ml}/10 \text{ ml} = 0.01 \text{ or } 1/100$$

The **Dilution Factor (DF)** is the reciprocal of the dilution, or **1/dilution**.

The **Total Dilution** is the product of all dilutions:

$$D_{total} = D_1 \times D_2 \times D_3 \dots$$

What does all this mean for spread plates and pour plates? The initial concentration, C_1 (that is, the concentration of cells in the undiluted culture) is equal to C_2 (the number of colonies/volume plated, or P_2/V_2) x DF, or:

$$C_1 = C_2 DF = (\frac{P_2}{V_2}) DF$$

where

 C_1 = initial concentration (or concentration in the stock culture)

 $C_2 = CFU/ml$ (or the concentration of CFU on the plate)

 $P_2 = CFU$ on the plate

V₂= Volume plated (ml)

DF = Dilution Factor

Dilution Problems

- 1. If there are 2×10^{12} bacteria in 200 ml, what is the cell concentration?
- 2. If you have 175 ml of a crystal violet solution of 5 µg/ml, how much crystal violet do you have?
- 3. If you add 2 ml of a bacterial suspension with a concentration of 2 x 10⁹ cells/ml to 13 ml of water, what is the new cell concentration?
- 4. Diagram a scheme to make a 1:10,000 dilution
- 5. Diagram a scheme to make a 1:7.2 dilution
- 6. You perform the following series of dilutions: 1:10, 1:3, 1:2. What is your total dilution? What is the DF?
- 7. After diluting your culture 1:5000, you have a cell concentration of 2.3×10^2 cells/ml. What was the initial concentration?
- 8. After diluting your culture 1:2500, you plate 0.1 ml and get 154 colonies. What was the initial concentration?

Dilution Problems (solutions)

1. If there are 2×10^{12} bacteria in 200 ml, what is the cell concentration?

$$C = \frac{P}{V} = \frac{(2 \times 10^{12} \text{ cells})}{(200ml)} = 1 \times 10^{10} \text{ cells /ml}$$

2. If you have 175 ml of a crystal violet solution of 5 μ g/ml, how much crystal violet do you have?

$$P = CV = (5 \mu g /ml)(175ml) = 875 \mu g crystal violet$$

3. If you add 2 ml of a bacterial suspension with a concentration of 2 x 10⁹ cells/ml to 13 ml of water, what is the new cell concentration?

$$C_1V_1 = C_2V_2$$

$$(2 \times 10^9 \text{ cells /ml})(2 \text{ ml}) = C_2(15 \text{ ml})$$

$$C_2 = \frac{4 \times 10^9 \text{ cells}}{15 \text{ ml}} = 2.67 \times 10^8 \text{ cells /ml}$$

4. Diagram a scheme to make a 1:10,000 dilution

$$\frac{1 \ ml}{1 \ ml + 9 \ ml} \rightarrow \frac{1}{1 + 9} \rightarrow \frac{1}{1 + 9} \rightarrow \frac{1}{1 + 9}$$
Total Dilution:
$$\frac{1}{10} \qquad \frac{1}{100} \qquad \frac{1}{1000} \qquad \frac{1}{10000}$$

5. Diagram a scheme to make a 1:7.2 dilution

$$\frac{1 \ ml}{1 \ ml + 6.2 \ ml} = \frac{1}{7.2}$$

6. You perform the following series of dilutions: 1:10, 1:3, 1:2. What is your total dilution? What is the DF?

$$Dil_t = \frac{1}{10} \times \frac{1}{3} \times \frac{1}{2} = \frac{1}{60}$$

$$DF = \frac{1}{dilution} = \frac{1}{(\frac{1}{60})} = 60$$

7. After diluting your culture 1:5000, you have a cell concentration of 2.3×10^2 cells/ml. What was the initial concentration?

$$C_1 = C_2DF = (2.3 \times 10^2 \text{ cells/ml})(5000) = 1.15 \times 10^6 \text{ cells/ml}$$

8. After diluting your culture 1:2500, you plate 0.1 ml and get 154 colonies. What was the initial concentration?

$$C = (\frac{P}{V})DF = (\frac{154 \ CFU}{0.1 \ ml})(2500) = 3.85 \ x \ 10^6 \ CFU/ml$$

Dilution Problems (Set #2)

- 1. Diagram a scheme to make a 1:3500 dilution. (Just draw the dilutions, not the tubes.)
- 2. You dilute a culture of *Saccharomyces cerevisiae* 1:100,000 and plate 0.2 ml. After incubation, you count 60 colonies. What was the concentration of cells in the original culture?
- 3. A spread plate inoculated with 0.2 ml from a 10⁻⁸ dilution contained 90 colonies. Calculate the cell concentration of the original culture.
- 4. Your spread plates contain 23, 191, and 643 CFU for your 10⁻⁶, 10⁻⁷, and 10⁻⁸ dilutions, respectively. What was the concentration of cells in the original culture?
- 5. Starting with a stock bacterial culture with a **total** cell count of 6.2 x 10¹⁰ cells/ml, outline an experiment for performing a viable count using pour plates. Assume that at least 10% of the culture is still viable. Provide:
 - a. a dilution scheme for the stock culture, including the volumes of cells and diluent used
 - b. the dilutions used and the volume plated
 - c. the rationale for using the dilutions selected
- 6. You decide to assay a stock solution of bacteria that is supposed to have 3.7×10^{8} cells/ml. In order to perform the assay, show:
 - a. how you would dilute the stock (volumes of stock and diluent)
 - b. the dilutions you would plate, as well as the volume plated
 - c. the number of cfu expected

Dilution Problems (Set #2) SOLUTIONS

1. Diagram a scheme to make a 1:3500 dilution. (Just draw the dilutions, not the tubes.)

dilution 1/10 1/10 1/10 1/3.5 1/10 = 1 ml culture + 9 ml diluent 1/3.5 = 1 ml culture + 2.5 ml diluent

total dilution 10^{-1} 10^{-2} 10^{-3} 2.86 x 10^{-4} (or 1/3500)

DF 10 100 1000 3500

2. You dilute a culture of *Saccharomyces cerevisiae* 1:100,000 and plate 0.2 ml. After incubation, you count 60 colonies. What was the concentration of cells in the original culture?

$$C_1 = C_2 \times DF = (P/V) \times DF = (60 \text{ cfu}/0.2 \text{ ml}) \times 100,000 = 3.0 \times 10^7 \text{ cfu/ml}$$

3. A spread plate inoculated with 0.1 ml from a 10^{-8} dilution contained 90 colonies. Calculate the cell concentration of the original culture.

$$C_1 = C_2 \times DF = (P/V) \times DF = (90 \text{ cfu/0.1 ml}) \times 10^8 = 9.0 \times 10^{10} \text{ cfu/ml}$$

4. Your spread plates contain 23, 191, and 643 CFU for your 10⁻⁶, 10⁻⁷, and 10⁻⁸ dilutions, respectively. What was the concentration of cells in the original culture?

First, select the proper plate to use. We want between 30 and 300 colonies/plate, so that would be the one with 191 colonies, which was inoculated from the 10⁻⁷ dilution. Since the volume plated wasn't given, we'll assume it was 0.1 ml, since that's most common with spread plates.

$$C_1 = C_2 \times DF = (P/V) \times DF = (191 \text{ cfu/0.1 ml}) \times 10^7 = 1.9 \times 10^{10} \text{ cfu/ml}$$

- 5. Starting with a stock bacterial culture with a **total** cell count of 6.2×10^{10} cells/ml, outline an experiment for performing a viable count using pour plates. Assume that at least 10% of the culture is still viable. Provide:
 - a. a dilution scheme for the stock culture, including the volumes of cells and diluent used
 - b. the dilutions used and the volume plated
 - c. the rationale for using the dilutions selected

We want 30-300 colonies/plate. With the total cell concentration given, 62 colonies is a good target number. If all of the cells are alive, the dilution that gives us 62 cfu in a 0.1 ml inoculum will have a concentration of 620 cfu/ml, since C = P/V = (62 cfu/0.1 ml).

 C_1 (stock culture concentration) = C_2 (concentration in the dilution plated) x DF

DF =
$$C_1/C_2$$
 = $(6.2 \times 10^{10} \text{ cfu/ml})/(620 \text{ cfu/ml})$ = 1×10^8

Total dilution = $1/DF = 10^{-8}$

- a. A dilution scheme to give a total dilution of 10^{-8} could be 1:100 \rightarrow 1:100 \rightarrow 1:10.
- b. Plate 0.1 ml (as used in the calculation above) from the 10⁻⁸ and 10⁻⁷ dilutions.
- c. We use two dilutions as inocula because the cells may only be 10% viable. That is, if the cells are all alive, the 0.1 ml from the 10-8 dilution will contain 62 cfu. But, if it's only 10% viable, it will contain only 1/10 as many cfu, or only about 6 cfu. If that's the case, we protect ourselves (since we won't know how many cells are viable for a couple of days) by plating from the next lower dilution. Since that dilution is 10x more concentrated (or less dilute) it will have 10x more viable cells, or 62 cfu/0.1 ml.

- 6. You decide to assay a stock solution of bacteria that is supposed to have 3.7×10^8 cells/ml. In order to perform the assay, show:
 - a. how you would dilute the stock (volumes of stock and diluent)
 - b. the dilutions you would plate, as well as the volume plated
 - c. the number of cfu expected

We want 30-300 colonies/plate. We could try for 37, but that's too close to 30, so we'll double it and aim for 74 colonies/plate. We could adjust for that in our dilution scheme, but it's easier to just double our inoculum to the maximum permitted volume of 0.2 ml. That means $C_2 = (74 \text{ cfu/0.2 ml}) = 370 \text{ cfu/ml}$.

DF =
$$C_1/C_2$$
 = (3.7 x 10⁸ cfu/ml)/(3.7 x 10² cfu/ml) = 1 x 10⁶; Total dilution = 10⁻⁶

a. For a total dilution of 10^{-6} we could do serial 1/10 dilutions = 1/10 \rightarrow 1/10 \rightarrow 1/10 \rightarrow 1/10 \rightarrow 1/10

We do 1/10 dilutions (instead of some 1/100) because we have no idea how many cells are viable and will need to plate from lots of dilutions in case viablility is very low. (Remember, with 3.7×10^8 cells/ml if even 0.0001% are alive, that's still 3.7×10^2 cells/ml.)

- b. As indicated above, since we don't know the concentration of viable cells, we plate from a number of dilutions in order to cover a wide range of possibilities. I would probably use every dilution, just to be safe. That would cover me from viable cell concentrations of 370 cfu/ml up through 3.7 x 10⁸ cfu/ml.
- c. Since I would be plating 0.2 ml, I would expect 74 cfu from the 10⁻⁶ dilution if all cells were viable. As viability dropped, each lower dilution would compensate for decreases of 10%. As I stated above, by plating 0.2 ml from each dilution, I would be covered for 100% viability, down to 0.0001%.