



**dpl<sup>®</sup>**

Neck Pain Relief

Light Therapy System

## INCLUDED

- **dpl**<sup>®</sup> Neck Pain Relief
- User Guide
- Power Supply
- Power Cable

LED Lights



Neck Pillow

# WELCOME

Welcome to LED Technologies, Inc and our **dpl**<sup>®</sup> Neck Pain Relief System. We are thrilled to announce this product as an in-home device to be used at your convenience.

The **dpl**<sup>®</sup> Neck Pain Relief System uses light emitting diodes (LEDs) in the infrared and red spectrums, which are medically proven to relax muscles and relieve muscle spasms, pain and stiffness.

Please read this manual before operating the **dpl**<sup>®</sup> Neck Pain Relief System.

## CAUTIONS

1. The **dpl**<sup>®</sup> Neck Pain Relief System is an electrical device. To avoid electrical shock:
  - Do not wash electrical parts with water.
  - Do not touch the plug or device with wet hands.
  - Do not touch unit to wet skin.
  - Do not use if cord or plug appears to be damaged.
  - Do not use the unit if it has been damaged, torn, or has come in contact with water.
  - Do not pull, lift, or carry unit by the power cord.

1. Do not use the cord as a handle, close a door on the cord, or pull the cord against sharp edges or corners.
2. Do not use near heated surfaces.
3. Do not operate the unit in areas where flammable or combustible products or vapors may be present. Explosion or fire may occur.
4. Do not allow children to use – this is not a toy. Do not leave small children unattended close to unit.
5. Not to be used by disabled, sleeping or unconscious persons.
6. Do not use the unit on an open wound, as heat may increase circulation and cause bleeding.
7. Do not use with liniments, salves or ointments that contain heat-producing ingredients as it may cause skin discomfort.

**Consult your doctor prior to using, especially under these circumstances:**

- You are pregnant, are taking drugs known to cause sensitivity of skin to light, have any suspicious or cancerous lesions, and/or recently had a steroid injection or use topical steroids.
- Read all instructions and safety precautions before using.
- If you experience discomfort or have a concern about its use, stop use immediately and contact your doctor.
- Use only the supplied power supply.

- Unplug from outlet when not in use to avoid tripping hazards.
- Do not unplug by pulling the cord. Grasp the plug itself, not the cord, when unplugging from an outlet.
- Do not exceed 40 minutes of use in one therapy.

## OPERATING YOUR NECK PAIN RELIEF SYSTEM

- Plug device USB cord into power supply or other universal USB power port.
- Place device on neck.
- Turn on device with on/off button on power cord.
- Enjoy the therapeutic lights.
- Treatment time is 14-25 minutes.
- Turn off device and unplug the adapter after treatment.

## TREATMENT

- There is no timer. Treatment time is managed by user.
- Treatment time is 14-25 minutes.
- Use daily for maximum results.

## HOW TO CLEAN DEVICE

- Wipe clean with damp cloth or alcohol wipe.
- Cover can be removed and washed.
- To wash cover, disconnect power cord. Unzip cover and gently remove to ensure no damage the LED lights.
- Install cover gently after washing, ensure all lights and the power cord are exposed. Zip up zipper.

## CUSTOMER SERVICE

Need further assistance? You can reach LED Technologies, Inc's Customer Service by calling (800) 337-9565, Monday-Friday, 7:30am-4:00pm, Mountain Standard Time, or by emailing [customerservice@ledtechnologies.com](mailto:customerservice@ledtechnologies.com).

## CONNECT WITH US

When you register your **dpl**<sup>®</sup> Neck Pain Relief System at [revivelighttherapy.com/warranty](http://revivelighttherapy.com/warranty), you will be entered into our system to be eligible for exclusive offers, VIP discounts, new product previews, and the opportunity to earn free products through new product trials.

## REPLACEMENT PARTS

The following lost or damaged parts can be purchased as replacement items directly from LED Technologies, Inc by calling Customer Service at (800) 337-9565, or by emailing [customerservice@ledtechnologies.com](mailto:customerservice@ledtechnologies.com).

- Power Supply
- Power Cord

## DISCLAIMER

This product is not intended to diagnose, cure or prevent disease. LED Technologies, Inc makes no claims, representations or warranties regarding the ability of this product to cure any physical, oral or mental condition. A qualified health professional should be consulted with regard to any condition requiring medical attention.

## WARRANTY

LED Technologies, Inc warrants that this product will be free from defects in materials and workmanship for 1 year from the date of purchase by the original owner. This warranty covers the original purchaser only and is not assignable or transferable to another party. Register your warranty online at [revivelighttherapy.com/warranty](http://revivelighttherapy.com/warranty)



6000 Greenwood Plaza Blvd Ste 110  
Greenwood Village, CO 80111  
800.337.9565

Our products are continually  
improving. The actual product may  
differ from the product shown.