

**LC  
S**



National Drawdown Framework Agreement for the Provision of

# **AIDS & APPLIANCES**

Therapeutic Aids - HSE Ref: 9231



Simple  
Stuff  
Works

From our family to yours

[lcs seating.ie](http://lcs seating.ie)



## About Simple Stuff Works®

We are a family run organisation with a passion for supporting people to protect body shape. We provide high quality education and therapeutic positioning equipment with an aim to eradicate the preventable problems associated with body shape distortion. We are driven by our core family values of compassion, trust and dignity. Our advisors and colleagues are all highly experienced, having written some of the leading research into Postural Care and its effects.

Our equipment is suitable for those with relatively uncomplicated or the most complex of healthcare needs. We pride ourselves on the fact that our equipment is silent, hygienic and simple to use whilst our training programmes will allow you to implement this powerful form of therapy as safely and as humanely as possible.

## CONTENTS

Paediatric Night Positioning System (LOT NO. 6.1)	04
Paediatric Day Positioning System (LOT NO. 6.2)	06
Adult Night Positioning System (LOT NO. 6.3)	08
Adult Day Positioning System (LOT NO. 6.4)	10
Additional products	12
Additional products	14

# Paediatric Night Positioning System



## Top-to-Toe Stabilising Mesh

The base of the award-winning sleep system starts with two layers of non-slip mesh. This can be either one layer of fitted mesh with a flat top piece, or two flat pieces if the individual is maybe sharing a bed. The flexibility of this base layer makes it one of the most dynamic parts of our system.

*Single Fitted Mesh:*  
 Bottom piece - 1830mm x 2500mm  
 Top Piece - 1830mm x 1000mm  
*Cot Fitted Mesh:*  
 Bottom piece - 1830mm x 1300mm  
 Top piece - 1500mm x 700mm



## Mattress Topper

The Temperate Topper is a soft mattress topper which is used underneath the sheet on top of the layers of Stabilising Mesh. The topper wicks moisture away from the skin and allows air to circulate around an individual, helping to maintain a dry, temperate environment for them to rest in. The Temperate Topper is non-dynamic, neither cooling nor warming the individual.

*Available in single and cot size*



## Fitted Sheets

Our sheets are made from a beautiful material which is hard wearing but lovely and soft. They are extra big and stretchy so that supports can be placed around an individual and the sheet will still conform to their body.

*Cot bed: 700mm x 1400mm*  
*Single bed: 900mm x 1900mm*  
*Double bed: 1350mm x 1900mm*



## Padded Hip Supports

The Padded Hip Supports providing firmer support for an individual. They are L shaped with soft, washable, removable pads which cover them. These are then positioned in between the two layers of mesh which hold them in place securely.

*Small - 100mm x 100mm x 100mm*  
*Medium - 140mm x 140mm x 140mm*  
*Large - 200mm x 200mm x 200mm*



## Padded Trunk Supports

The padded lateral supports are perfect for providing firmer support for an individual using our therapeutic positioning equipment. They are L shaped with soft, washable, removable pads which cover them. These are then positioned in between the two layers of mesh which hold them in place securely. They can also be used to support other items such as legs supports or pillows.

*Small - 100mm x 100mm x 100mm*  
*Medium - 140mm x 140mm x 140mm*  
*Large - 200mm x 200mm x 200mm*



## Size 1 Supine Stabilizer

The Supine Stabilizer is a leg support piece which is designed to help an individual to keep a neutral, supine position for long periods of time. It is worn above the knees, across the thighs, and it stabilizes the pelvis and the femura.

*Size 1 = 450mm width x 250mm height x 130mm depth*  
*Size 2 = 450mm width x 300mm height x 150mm depth*  
*Size 3 = 500mm width x 300mm height x 200mm depth*  
*Size 4 = 600mm width x 350mm height x 225mm depth*  
*Size 5 = 750mm width x 350mm height x 225mm depth*



Seirbhís Sláinte  
 Níos Fearr  
 & Forbairt

Building a  
 Better Health  
 Service

**LOT NO. 6.1**

# Paediatric Day Positioning System



## Bed in a Bag

A Bed in a Bag comprises of a zip closing carry bag, 10cm thick fibre mattress, top-to-toe stabilising mesh, mattress topper and a fitted sheet. There is also the option to have a Treat Eeiz topper. The bed in a bag can be easily folded and stored away when not in use. It can be used as a mobile mattress. A bed in a bag is very useful for sensory rooms.

COT - 120cm x 60cm  
Cot Bed - 140cm x 70cm  
Mat - 190cm x 91cm



## Padded Hip Supports

The padded hip supports providing firmer support for an individual. They are L shaped with soft, washable, removable pads which cover them. These are then positioned in between the two layers of mesh which hold them in place securely.

Small - 100mm x 100mm x 100mm  
Medium - 140mm x 140mm x 140mm  
Large - 200mm x 200mm x 200mm



## Size 1 Supine Stabilizer

The Supine Stabilizer is a leg support piece which is designed to help an individual to keep a neutral, supine position for long periods of time. It is worn above the knees, across the thighs, and it stabilises the pelvis and the femura.

Size 1 = 450mm width x 250mm height x 130mm depth  
Size 2 = 450mm width x 300mm height x 150mm depth  
Size 3 = 500mm width x 300mm height x 200mm depth  
Size 4 = 600mm width x 350mm height x 225mm depth  
Size 5 = 750mm width x 350mm height x 225mm depth



## Padded Trunk Supports

The padded lateral supports are perfect for providing firmer support for an individual using our therapeutic positioning equipment. They are L shaped with soft, washable, removable pads which cover them. These are then positioned in between the two layers of mesh which hold them in place securely. They can also be used to support other items such as legs supports or pillows.

Small - 100mm x 100mm x 100mm  
Medium - 140mm x 140mm x 140mm  
Large - 200mm x 200mm x 200mm



LOT NO. 6.2



# Adult Night Positioning System



## Top-to-Toe Stabilising Mesh

The base of the award-winning sleep system starts with two layers of non-slip mesh. This can be either one layer of fitted mesh with a flat top piece, or two flat pieces if the individual is maybe sharing a bed. The flexibility of this base layer makes it one of the most dynamic parts of our system.

*Single Fitted Mesh:*  
 Bottom piece - 1830mm x 2500mm  
 Top Piece - 1830mm x 1000mm  
*Cot Fitted Mesh:*  
 Bottom piece - 1830mm x 1300mm  
 Top piece - 1500mm x 700mm



## Padded Hip Supports

The padded hip supports providing firmer support for an individual. They are L shaped with soft, washable, removable pads which cover them. These are then positioned in between the two layers of mesh which hold them in place securely.

*Small - 100mm x 100mm x 100mm*  
*Medium - 140mm x 140mm x 140mm*  
*Large - 200mm x 200mm x 200mm*



## Size 3 Supine Stabilizer

The Supine Stabilizer is a leg support piece which is designed to help an individual to keep a neutral, supine position for long periods of time. It is worn above the knees, across the thighs, and it stabilizes the pelvis and the femora.

*Size 1 = 450mm width x 250mm height x 130mm depth*  
*Size 2 = 450mm width x 300mm height x 150mm depth*  
*Size 3 = 500mm width x 300mm height x 200mm depth*  
*Size 4 = 600mm width x 350mm height x 225mm depth*  
*Size 5 = 750mm width x 350mm height x 225mm depth*



## Padded Trunk Supports

The padded lateral supports are perfect for providing firmer support for an individual using our therapeutic positioning equipment. They are L shaped with soft, washable, removable pads which cover them. These are then positioned in between the two layers of mesh which hold them in place securely. They can also be used to support other items such as legs supports or pillows

*Small - 100mm x 100mm x 100mm*  
*Medium - 140mm x 140mm x 140mm*  
*Large - 200mm x 200mm x 200mm*



LOT NO. 6.3

# Adult Day Positioning System



## Bed in a Bag

A Bed in a Bag comprises of a zip closing carry bag, 10cm thick fibre mattress, top-to-toe stabilising mesh, mattress topper and a fitted sheet. There is also the option to have a Treat Eeiz topper. The bed in a bag can be easily folded and stored away when not in use. It can be used as a mobile mattress. A bed in a bag is very useful for sensory rooms.

COT - 120cm x 60cm  
Cot Bed - 140cm x 70cm  
Mat - 190cm x 91cm



## Padded Hip Supports

The Padded Hip Supports providing firmer support for an individual. They are L shaped with soft, washable, removable pads which cover them. These are then positioned in between the two layers of mesh which hold them in place securely.

Small - 100mm x 100mm x 100mm  
Medium - 140mm x 140mm x 140mm  
Large - 200mm x 200mm x 200mm



## Size 3 Supine Stabilizer

The Supine Stabilizer is a leg support piece which is designed to help an individual to keep a neutral, supine position for long periods of time. It is worn above the knees, across the thighs, and it stabilizes the pelvis and the femura.

Size 1 = 450mm width x 250mm height x 130mm depth  
Size 2 = 450mm width x 300mm height x 150mm depth  
Size 3 = 500mm width x 300mm height x 200mm depth  
Size 4 = 600mm width x 350mm height x 225mm depth  
Size 5 = 750mm width x 350mm height x 225mm depth



## Padded Trunk Supports

The padded lateral supports are perfect for providing firmer support for an individual using our therapeutic positioning equipment. They are L shaped with soft, washable, removable pads which cover them. These are then positioned in between the two layers of mesh which hold them in place securely. They can also be used to support other items such as legs supports or pillows.

Small - 100mm x 100mm x 100mm  
Medium - 140mm x 140mm x 140mm  
Large - 200mm x 200mm x 200mm



LOT NO. 6.4



### The Side Lying Leg Support

Our Side Lying Leg Support is designed to protect the hips and pelvis when a person is lying on their side. The leg channel creates a perfect place for the top leg to lie and the bottom leg wraps gently around the outside. When combined with lots of support at the back of the person, a lovely supportive side lying position can be achieved. This is perfect for individuals who need to lie on their side for safety reasons. The Side Lying Leg Support comes in three different sizes and a choice of three colours.

*Size 1 - 380mm length x 300mm depth with a leg groove of 120mm at widest point and 50mm at narrowest*

*Size 2 - 600mm length x 400mm depth with a leg groove of 120mm at widest point and 80mm at narrowest*

*Size 3 - 600mm length x 400mm depth with a leg groove of 140mm at widest point and 100mm at narrowest*



### Support Options

We supply a wide range of options for making supports within our therapeutic positioning system. We can offer wedges, rolls and padded brackets which can be used in any combination to make a comfortable position for the individual.



### Foot Supports

Our foot supports are soft supports which help to hold the feet in a neutral position for long periods of time. They are made from the same fibre as our Temperate Toppers so they help to keep the feet feeling fresh whilst they are protected from falling forwards or rotating. The foot supports come in three sizes and three colours.

*Size 1 - fits from baby up to adult size 1*

*Size 2 - fits up to adult size 5*

*Size 3 - fits up to size 9*



### Pillows

We have lots of different pillow options which can be used to support an individual. We offer regular pillows, horseshoe pillows, sausage pillows and neck support pillows. These all come in three different colours and the firmness of each pillow can be adapted depending on the needs and comfort of the individual.



Seirbhís Sláinte  
Níos Fearr  
& Forbairt | Building a  
Better Health  
Service

**LOT NO. 6.4**



### Treat Eeiz Pressure Care Relief

The Treat Eeiz full length overlay is a revolutionary mattress topper which can be used on any standard single or double bed without affecting the partners sleep. Incredibly, it has been clinically proven to help both prevent and heal stage 1, 2, 3 and 4 pressure ulcers by Leading Tissue Viability Nurses.

The Treat-Eeezi is far more comfortable than plastic or nylon air inflated products and there is no electric pump noise going on all night. It also cannot be punctured, it cannot fail in the middle of the night and it will not bottom out.

The Treat-Eeezi decreases sweating allowing the patient far superior sleep and protection even in a power failure with zero maintenance costs, zero running costs, no expert training required, no pumps to fail and student nurse friendly.

The Treat-Eeezi reduces shearing like no other product and the patient lies in a static system without repeated movement throughout the night. It can be washed at home to 80 degrees, bleached if stained and also tumbled dried at low heat in time for same evening use.

For patients who may not want to use electric air beds the Treat-Eeezi is a must. It can also be custom made for double beds or for bariatric patients up to 50 stone in weight.



### Multi-Plane Stabilizer

Our New Multi-Plane Stabilizer is a beautifully dynamic leg support piece which supports an individual's legs in a therapeutic position.

It has the ability to either internally or externally rotate the legs depending on what position would be beneficial for the individual, as well providing a nice level of abduction of the hips to keep them comfortable and in joint. It is also multi-sized, meaning it can grow with someone and does not need to be replaced as often.

There are two sizes of multi plain supine stabilizer size 1 and 2. Within these sizes there are different size options, offering two sizes within a size.

Colour options - cream, pink and blue  
Material options - a wipe down material is also available.







*“The LCS range of custom and combination seating accessories can be added to all our chairs, further tailoring it to the needs of the user.”*

*Please contact your local LCS clinical sales member for more details.*

.....

**Luke** (Donegal)

phone: **071 9163905** email: [luke@lcseating.ie](mailto:luke@lcseating.ie)

**Grace** (Connaught)

phone: **086 460 7954** email: [grace@lcseating.ie](mailto:grace@lcseating.ie) office: **071 9163905**

**Richard** (Leinster)

phone: **086 410 4762** email: [richard@lcseating.ie](mailto:richard@lcseating.ie) office: **071 9163905**

**Kevin** (Munster)

phone: **086 170 0499** email: [kevin@lcseating.ie](mailto:kevin@lcseating.ie) office: **071 9163905**

**David** (Ulster/NI North Leinster)

phone: **00353 86 836 6973** or **0044 77 926 94653** email: [david@lcseating.ie](mailto:david@lcseating.ie)  
office: **071 9163905**

.....

General queries: [info@lcseating.ie](mailto:info@lcseating.ie)

Appointments: [appointments@lcseating.ie](mailto:appointments@lcseating.ie)

**lcseating.ie**

Unit 1, Grange Business Park,  
Grange, Co. Sligo,  
Ireland, F91 K300

Email: [info@lcseating.ie](mailto:info@lcseating.ie)  
Tel: 071 9163905  
Fax: 071 9163918